

## **DNAgard® Saliva HT Frequently Asked Questions (FAQ's)**

- 1. What is DNAgard® Saliva HT?**
- 2. Can DNAgard® Saliva HT be used for diagnostics purposes?**
- 3. What features make DNAgard® Saliva HT donor friendly?**
- 4. What are the specifics of the tracking labels?**
- 5. How long is the saliva DNA sample stabilized upon mixing with the DNAgard Saliva HT solution?**
- 6. What volume of saliva is needed for DNAgard® Saliva HT?**

### **1. What is DNAgard® Saliva HT?**

DNAgard® Saliva HT is a saliva DNA collection device optimized for laboratory processing with automated instruments. Its pierce-able cap and triple redundant labeling make it ideal for laboratories performing high-throughput DNA analysis. Additionally, DNAgard® Saliva HT features a unique, user-friendly design, with a cap-integrated stabilizer and leak prevention features to improve donor compliance.

### **2. Can DNAgard® Saliva HT be used for diagnostics purposes?**

No. DNAgard® Saliva HT is For Research Use Only. Not for use in diagnostic procedures.

### **3. What features make DNAgard® Saliva HT donor friendly?**

The device features integrated stabilizer, which is released from the funnel lid into the saliva sample in the tube when donors screw on the lid. Additionally, the device features a wide funnel for cleaner saliva collection, integrated stabilizer to avoid user contact, a tethered cap to prevent loss, and a wide base to prevent spills.

### **4. What are the specifics of the tracking labels?**

The triple redundant labeling, which is CDC compliant and LIMS compatible includes a laser etched barcode on the side of the tube, a 2D barcode on the bottom of the tube, and a readable label on the side of the tube.

### **5. How long does DNAgard® Saliva HT stabilize saliva DNA?**

1 year.

### **6. What volume of saliva needs to be collected for DNAgard® Saliva HT?**

2 mL.

#### **BIOMATRICA, INC.**

5627 Oberlin Drive, Suite 120,  
San Diego, California 92121 USA

Tel: 1-866-379-6879 • 858-550-0308 • Fax: 858-678-0597

Email: [info@biomatrica.com](mailto:info@biomatrica.com) • Web: [www.biomatrica.com](http://www.biomatrica.com)